

vernachair

The versatile 3-in-1
hygiene chair system

Commode
Shower Chair
Sani-chair

Product Code: 640ZA001



 Vernacare



From the bedside to the privacy of the bathroom, the versatile 3-in-1 hygiene chair system is designed for safe and comfortable patient toileting and washing.

Its modern smooth seamless design reduces the harbouring of micro-organisms and ensures greater patient comfort. Compatible with disposable bedpans and commode pans, the chairs are integral to the Vernacare single-use system, used in over 94% of UK hospitals.

Developed in collaboration with clinicians from across the healthcare industry and Independently tested (compliant to BS 4751).

VernaChair Specification

Height	1081mm
Width	607mm
Depth	845mm
Width (between arms)	480mm
Seat height	570mm
Maximum patient weight	159kg/25st

VC300

For more information:

t: +44 (0)1204 529494 e: nathan.shields@vernagroup.com w: vernacare.com
Vernacare Ltd | Folds Road | Bolton | Lancashire | BL1 2TX | UK

Vernacare