

VernaChair

VernaChair can be used as a commode chair, a showerchair or a sanichair.

VernaChair

VernaChair
SHOWER

VernaChair can be quickly switched from a **Commode** to a **Shower** chair by changing the seat!



- ✓ **Only 5 removable parts**
for the VernaChair Commode
- ✓ **Only 2 removable parts**
for the VernaChair Shower
- ✓ **Easy to assemble**
no tools required
- ✓ **Covered castors**
help prevent debris collection
- ✓ **Stainless steel frame**
withstands modern cleaning methods
- ✓ **Independently tested**
and compliant to BS 4751

**Ideal
for
infection
prevention
and control**


Vernacare

Developed in collaboration with clinicians from across the healthcare industry.


From the bedside to the privacy of the bathroom, VernaChair is designed for safe and comfortable patient toileting and washing. Its modern smooth seamless design reduces the harbouring of micro-organisms and ensures greater patient comfort.

The **VernaChair commode** dismantles into just **five parts** for easy cleaning and reassembly without the use of tools. The chair's sturdy stainless steel frame and covered castors are backed by a 3 year warranty.

Compatible with disposable bedpans and commode pans, the VernaChair commode is integral to the Vernacare single-use system, used in over 94% of UK hospitals.



Rear loading pan for easy removal



Removable footrest for easy cleaning



Hinged moveable armrests for confident patient transfer



Large easy-locking castors make the chair easy to manoeuvre and to secure

The **VernaChair shower** is designed for comfortable and easy patient washing and showering. Its lightweight and clip-in seat is chemically resistant and allows quick and easy draining. With only **two removable parts**, the VernaChair shower may be disassembled quickly and easily for thorough cleaning, reducing the harbouring of micro-organisms.

We recommend not to use the footrest with the VernaChair shower as it may be slippery.



Clip-in seat

Cleaner by Design.

VernaChair Specifications	
Full Height	1070mm
Width	600mm
Depth (front to back)	814mm
Width between Arms	475mm
Seat height from floor (Commode)	562.5mm
Height from floor to seat base (Shower)	567mm
Height from floor to seat top (Shower)	600mm
Maximum patient weight	159kgs / 25 st

NPC Code: VernaChair Commode GTB1563

Product Codes: VernaChair Commode 600ZA001
VernaChair Shower Chair 630ZA001
Shower Chair Seat 631ZA001

Spare parts are also available.

Vernacare

1 Western Avenue,
Matrix Park,
Chorley PR7 7NB
T: +44 (0)1204 529494
E: info@vernacare.com
www.vernacare.com

VC004 | REV00 | May 2018