
ELECTRONIC FLUID MEASURING SCALE



ELECTRONIC FLUID MEASURING SCALE

- ✦ Hygienic
- ✦ Easy to use
- ✦ Maximum Weight: 5kg/11lb
- ✦ 12 Month Warranty (excludes misuse)



Code	Description	Average Weight
101AA100	Bedpan	47.9 +/- 5 g
102AA100	Male Urinal (800ml)	37.5 +/- 5g
103AA120	Male Urinal (900ml)	37.0 +/- 6g
104AA200	General Purpose bowl	21.0 +/- 3g
105AA260	Kidney Bowl	19.2 +/- 3g
106AA300	CS Receiver	24.8 +/- 6g
108AA024	VernaFem	58.0 +/- 5g
109AA100	Midi Slipper Pan	46.5 +/- 5g
110AA100	Jug	27.5 +/- 3g
111AA200	Medium Bowl	39.0 +/- 5g
112AA100	Large Bowl	48.0 +/- 5g
114AA100	Maxi Slipper Pan	51.0 +/- 5g
116AA200	Multi Cup	8.1 +/- 1g
118AA100	Washbowl	58.5 +/- 3g
121AA075	Bedpan (new)	50.0 +/- 5g
122AA050	Disposable Support	53.0 +/- 5g
123AA150	Commode Pan	39.0 +/- 5g

This chart is for guidance only. Due to natural composition of the material used in our products, all quoted weights are approximate and subject to fluctuation. As with all products, use and determination remains the responsibility of the user. Note when weighing 1ml = 1gm. When using VernaGel, please remember that each sachet weights 7g.

Product code: 599ZS001