

Midi Slipper Pan
109AA100
NPC Code FSE 014
Quantity 100 per Carton



KM 650633



Tested to Comply with
 PAS 029:1999

- SINGLE USE
- MOULDED FIBRE PULP
- MACERATABLE
- MADE FROM RECYCLED PRODUCT

- ☐ Ideal for general toileting of patients with limited mobility, in bed.
- ☐ Wedge-shape design slips underneath patient without lifting or rolling.
- ☐ Disposable single use product to reduce risk of cross infection
- ☐ Moulded fibre products made from 100% recycled paper.
- ☐ Fully disposable in Vernacare Disposal Units

Regulatory Compliance

- ☐ Compliant with: Medical Device Directive 93/42/EC Class 1
- ☐ Product Specification Assessment PAS 29:1999
- ☐ Disposability: By Maceration
- ☐ ISO 9001:2015 / ISO 14001:2015

Description:

- ☐ Capacity: 1.3 litre (max)
- ☐ Material: 100% Paper Pulp
- ☐ Dimensions: 330 x 290 x 85mm
- ☐ Weight: 46.5g (+/- 5g)
- ☐ Colour: Grey

Fluid Retention Testing

- ☐ Tested in Accordance with PAS 29:1999
- ☐ Guaranteed to hold fluid for 4 Hours
- ✓ Resistance Temperature: 35 C (+/- 3)



Outer Case Dimensions 350 x 605 x 295mm



OUTER BAR CODE



Due to the nature of the manufacturing process the products may vary in appearance. This will not affect the performance of the product.