

**Large Bowl
112AA100
NPC Code FSL 118
Quantity 100 per Carton**



KM 650633



Tested to Comply with
PAS29:1999

- SINGLE USE
- MOULDED FIBRE PULP
- MACERATABLE
- MADE FROM RECYCLED PRODUCT

- ☐ General purpose bowl.
- ☐ Ideal for handling larger volumes of liquid.
- ☐ Disposable single use product to reduce risk of cross infection
- ☐ Moulded fibre products made from 100% recycled paper.
- ☐ Fully disposable in Vernacare Disposal Units

Regulatory Compliance

- ☐ Compliant with: Medical Device Directive 93/42/EC Class 1
- ☐ Product Specification Assessment PAS 29:1999
- ☐ Disposability: By Maceration
- ☐ ISO 9001:2015 / ISO 14001:2015

Description:


- ☐ Capacity: 3.0 litre (max)
- ☐ Material: 100% Paper Pulp
- ☐ Dimensions: 300 x 100mm
- ☐ Weight: 48.0g (+/- 5g)
- ☐ Colour: Grey

Fluid Retention Testing

- ☐ Tested in Accordance with PAS 29:1999
- ☐ Guaranteed to hold fluid for 4 Hours
- ☐ Resistance Temperature: 35°C (+/- 3)



Outer Case Dimensions 300 x 580 x 300mm



(01)00800239011125

OUTER BAR CODE



Due to the nature of the manufacturing process the products may vary in appearance. This will not affect the performance of the product.