

# AZOWRAP™

## COMPLETE & COMPLETE HEAVY DUTY

The **AZOWRAP™** Complete and Complete Heavy Duty range offers the Complete wrapping solution for the busy healthcare SSD. Sheets of **AZOWRAP™** Non-woven and **AZOWRAP™** SMS material, bonded along two sides provide a quick, efficient, and aseptic single-step wrapping system, minimising overall processing time in the SSD.

### PRODUCT DESCRIPTION

- Bonded sheets of Non-woven and SMS material, latex free
- 57 gsm Inner Green non-woven layer, 55 gsm Outer Blue SMS layer (AZOWRAP™ Complete)
- 66 gsm Inner Green non-woven layer, 55 gsm Outer Blue SMS layer (AZOWRAP™ Complete Heavy Duty)
- Suitable for Steam (134°C/23 mins), Ethylene oxide (54°C/90 mins) and Plasma sterilisation processes
- 5 year Shelf life
- ISO11607-1 and EN868-2 compliant
- Single use only



### AZOWRAP™ COMPLETE

Product Code	Dimensions (mm)	Colour	Amount/Case
92275	900x900	Green/Blue	75
92276	1000x1000	Green/Blue	75
92277	1200x1200	Green/Blue	38
92278	1300x1500	Green/Blue	25
92279	1370x1370	Green/Blue	25
92280	1500x1800	Green/Blue	18

### AZOWRAP™ COMPLETE HEAVY DUTY

Product Code	Dimensions (mm)	Colour	Amount/Case
92281	900x900	Green/Blue	75
92282	1000x1000	Green/Blue	75
92283	1200x1200	Green/Blue	38
92284	1300x1500	Green/Blue	25
92285	1370x1370	Green/Blue	25
92286	1500x1800	Green/Blue	18

**AZOWRAP™**  
COMPLETE & COMPLETE HEAVY DUTY

**MATERIAL DATA** (SMS 55gsm single sheet)

Property	Units	Result (Typical)
Substance	g/m <sup>2</sup>	55
High Porosity	l/min/dm <sup>2</sup>	200
Tensile Strength MD	kN/m	2.5
Tensile Strength CD	kN/m	1.3
Wet Tensile Strength MD	kN/m	2.5
Wet Tensile Strength CD	kN/m	1.2
Burst Strength	kPa	245
Tearing Strength MD	mN	3500
Tearing Strength CD	mN	7000
Stretch MD	%	40
Hydrostatic test	mbar	30

**MATERIAL DATA** (Non-woven 57gsm single sheet)

Property	Units	Result (Typical)
Substance	g/m <sup>2</sup>	57
High Porosity	l/min/dm <sup>2</sup>	14
Tensile Strength MD	kN/m	2.5
Tensile Strength CD	kN/m	1.0
Wet Tensile Strength MD	kN/m	1.4
Wet Tensile Strength CD	kN/m	0.7
Burst Strength	kPa	190
Tearing Strength MD	mN	1100
Tearing Strength CD	mN	1550
Stretch MD	%	10
Hydrostatic test	mbar	30

**MATERIAL DATA** (Non-woven 66gsm single sheet)

Property	Units	Result (Typical)
Substance	g/m <sup>2</sup>	66
High Porosity	l/min/dm <sup>2</sup>	18
Tensile Strength MD	kN/m	2.5
Tensile Strength CD	kN/m	1.1
Wet Tensile Strength MD	kN/m	2.1
Wet Tensile Strength CD	kN/m	0.9
Burst Strength	kPa	240
Tearing Strength MD	mN	1350
Tearing Strength CD	mN	2050
Stretch MD	%	12
Hydrostatic test	mbar	29

**Also Available in the AZOWRAP™ Range**

- AZOWRAP™ Tray liner
- AZOWRAP™ Crepe
- AZOWRAP™ Non woven and Non woven Interleaved
- AZOWRAP™ SMS and SMS Interleaved
- AZOWRAP™ Combo Interleaved
- AZOWRAP™ Complete and Complete Heavy Duty**
- AZOWRAP™ Complete Plus And Complete Plus Heavy Duty
- AZOWRAP™ Transport
- AZOWRAP™ Transport Plus