

AZOWRAP™

TRAY LINER

AZOWRAP™ Tray liner can absorb up to 250% of its dry weight in water, helping to remove and dissipate residual moisture during the drying phase, and prevent wet packs.

AZOWRAP™ Tray liner can be used either inside, or outside the basket as an additional layer of protection.

PRODUCT DESCRIPTION

- Single layer of wet-laid cellulosic fibre material, latex free
- 70 gsm
- Suitable for use with Steam (134°C/23 mins) sterilisation processes
- 5 year shelf life
- ISO11607-1 (excluding Chap 5.2) compliant
- Single use only



Product Code	Dimensions (mm)	Colour	Amount/ Case
92218	350x250	White	1000
92219	300x500	White	500
92220	400x500	White	500
92221	560x480	White	500

MATERIAL DATA (Tray liner 60gsm, Single sheet)

Property	Units	Result (Typical)
Substance	g/m ²	70
High porosity	l/min/dm ²	15
Tensile Strength MD	kN/m	2.7
Tensile Strength CD	kN/m	1.4
Wet Tensile Strength MD	kN/m	0.9
Wet Tensile Strength CD	kN/m	0.5
Burst Strength	kPa	200
Tearing Strength MD	mN	1300
Tearing Strength CD	mN	1700
Stretch MD	%	9.5
Water Absorbency	%	250

AZOWRAP™

AZOWRAP™

TRAY LINER

Also Available in the AZOWRAP™ Range

AZOWRAP™ Tray liner

AZOWRAP™ Crepe

AZOWRAP™ Non woven and Non woven Interleaved

AZOWRAP™ SMS and SMS Interleaved

AZOWRAP™ Combo Interleaved

AZOWRAP™ Complete and Complete Heavy Duty

AZOWRAP™ Complete Plus And Complete Plus Heavy Duty

AZOWRAP™ Transport

AZOWRAP™ Transport Plus